

SAMPLING PROCEDURES- DRINKING WATER INORGANICS

Chain of Custody:

All samples must be submitted with an "Analysis Request Form / Chain of Custody".

Please fill in the form completely.

Keep the "PINK" copy of the chain of custody for your records

Return the chain of custody with your samples to the address below

CHEMTECH-FORD , Inc.
attn. Sample Receiving
6100 South Stratler (380 West)
Murray, UT 84107
Phone: 801-262-7299

General Instructions:

DO NOT rinse sample containers. Exercise caution while sampling.

Containers may contain an acid preservative necessary for testing.

Fill all sample containers to the neck of the bottle.

Wash your hands after sampling and be sure the lids are put back on tightly.

Please submit bottles on ice, keep ice in a Plastic Bag (4 degrees C).

Do Not Use Dry Ice. Bottles will freeze and break.

New Source Inorganics

Five Containers:

1- Unpreserved 1/2 Gallon

4- Cyanide Quart w/ NaOH (pink label)

2- Unpreserved Glass Amber Liter

5- Metals Pint w/ HNO3 (red label)

3- Nutrient Pint w/ H2SO4 (yellow label)

Inorganic & Metals

Three Containers:

1- Unpreserved Quart (If Nitrate or Nitrite is needed it will be taken out of this bottle).

2- Cyanide Quart w/ NaOH (pink label)

3- Metals Pint w/ HNO3 (red label)

Nitrate

One Container:

1- Unpreserved Pint

Radiological

One Container:

1- Radiological 1/2 Gallon w/ HNO3 (red label)

Sulfide

One Container:

1 - Sulfide Quart w/ Zinc Acetate and Sodium Hydroxide (gold label)

Holding Time

YOUR SAMPLES MUST BE RECEIVED AT CHEMTECH-FORD

WITHIN _____ HOURS OF SAMPLING.

Please do not submit samples with short holding times after 2:00 on Fridays,

Samples received at the lab after this time will need to be resampled.